

## Product Datasheet

### Coronavirus, Canine Antigen, Native LA235

<b>Product Name</b>	Coronavirus, Canine Antigen, Native
<b>Catalog Number</b>	LA235
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	≤ -70°C
<b>Note</b>	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
<b>Source</b>	A-72 cells/Virus strain 1-71
<b>Buffer</b>	0.05M Glycine, pH 9 with <0.1% CHAPS and <0.01% SDS
<b>Preservatives</b>	None