

## Product Datasheet

### Mycoplasma pneumoniae Grade 2 Antigen, Native LA240

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Mycoplasma pneumoniae Grade 2 Antigen, Native         |
| <b>Catalog Number</b>      | LA240   |
| <b>Tested Applications</b> | EIA   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Partially purified and concentrated by centrifugation |
| <b>Storage</b>             | $\leq -70^{\circ}\text{C}$                            |
| <b>Note</b>                | Sonicate immediately prior to usage                   |
| <b>Source</b>              | Strain FH   |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline                       |