

## Product Datasheet

### Human Immunodeficiency Virus 1 (HIV 1) Chimeric Antigen, Recombinant ≥80% LA249

<b>Product Name</b>	Human Immunodeficiency Virus 1 (HIV 1) Chimeric Antigen, Recombinant ≥80%
<b>Catalog Number</b>	LA249
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	2-8°C
<b>Purity</b>	≥80%
<b>Source</b>	
<b>Buffer</b>	Tris-HCl buffer, pH 8.5 with EDTA
<b>Specificity</b>	gp120 and gp41
<b>Matched Pairs</b>	LA038 (detection)
<b>Expression System</b>	E. coli