

Product Datasheet

Zopiclone-BSA Antigen LA269

Product Name	Zopiclone-BSA Antigen
Catalog Number	LA269
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	Distilled Water
Matched Pairs	HM310 (detection)
Preservatives	0.01% Sodium Azide