

Product Datasheet

EDDP-BSA Antigen LA282

Product Name EDDP-BSA Antigen

Catalog Number LA282

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Liquid; Antigen-BSA conjugate

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Matched Pairs HM381 (detection)

Preservatives 0.01% Sodium Azide