

Product Datasheet

Dehydroepiandrosterone sulfate (DHEA-S)-BSA Antigen LA285

Product Name	Dehydroepiandrosterone sulfate (DHEA-S)-BSA Antigen
Catalog Number	LA285
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM378, HM379 or HM380 (detection)
Preservatives	0.05% Proclin 950