

Product Datasheet

Ethyl glucuronide (EtG)-BSA Antigen LA299

Product Name	Ethyl glucuronide (EtG)-BSA Antigen
Catalog Number	LA299
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservatives	0.03% Sodium Azide