

Product Datasheet

Amphetamine-BSA Antigen LA311

Product Name	Amphetamine-BSA Antigen
Catalog Number	LA311
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	≥ 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.5
Matched Pairs	HM450 (detection)
Preservatives	0.01% Sodium Azide