

## Product Datasheet

### Neutrophil Gelatinase-associated Lipocalin (NGAL) Antigen, Recombinant >98% LA312

<b>Product Name</b>	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antigen, Recombinant >98%
<b>Catalog Number</b>	LA312
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (Full length with His-tag)
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>98% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.09% Sodium Azide