

## Product Datasheet

### Streptolysin O LA318

**Product Name** Streptolysin O

**Catalog Number** LA318

**Tested Applications** LETIA

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen

**Storage**  $\leq -20^{\circ}\text{C}$ . Avoid repeated freezing and thawing.

**Note** 10 mg minimum purchase quantity. Approximate antigenic activity of 4000 IU/ml and specific antigenic activity of 2000 IU/mg of protein.

**Purity** >90% by SDS-PAGE

**Buffer** Tris-HCl, pH 8.2

**Expression System** E. coli