

Product Datasheet

Streptolysin O LA318

Product Name Streptolysin O

Catalog Number LA318

Tested Applications LETIA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage $\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.

Note 10 mg minimum purchase quantity. Approximate antigenic activity of 4000 IU/ml and specific antigenic activity of 2000 IU/mg of protein.

Purity >90% by SDS-PAGE

Buffer Tris-HCl, pH 8.2

Expression System E. coli