

Product Datasheet

Norovirus genogroup I Antigen, Recombinant >95% LA320

Product Name	Norovirus genogroup I Antigen, Recombinant >95%
Catalog Number	LA320
Tested Applications	ELISA, Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (Full length VP1 peptide with His-tag at N-terminus)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline with 1M Urea and 50 mM Argine
Matched Pairs	HM637 or HM884
Expression System	E. coli
Preservatives	0.05% Sodium Azide