

Product Datasheet

Testosterone-3-BSA Antigen LA329

Product Name Testosterone-3-BSA Antigen

Catalog Number LA329

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Liquid; Antigen-BSA conjugate

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Purity ≥90%

Buffer 10 mM Phosphate Buffer, pH 7.2

Matched Pairs HM521, HM522 or HM523 (detection)

Preservatives 0.01% Sodium Azide