

Product Datasheet

Anti-h cTnl 9705 SPTN-1 100125

Product Name	Anti-h cTnl 9705 SPTN-1
Catalog Number	100125
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Alternative Names	troponin
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	1.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	In the N-terminal part of the cTnl molecule (residues 21-30)
Purity	≥ 95 %
Associated Products	Recombinant cTnl antigen 610102 Native cTnl antigen 550-11
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.4-7.5
Cross Reactivity	Does not recognize skeletal troponin I (sTnl)
Specificity	Antibody recognizes human cardiac troponin I (cTnl)
Matched Pairs	RC9750 is recommended as a very sensitive capture. For increased sensitivity, 2+1 or 2+2 antibody approach is recommended, for example capture: 9701/RC9701 & 9705, detection 9703 & 9707/RC9707.
Shelf Life	36 months