

Product Datasheet

Dexamethasone-BSA Antigen LA333

Product Name	Dexamethasone-BSA Antigen
Catalog Number	LA333
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM533, HM534, HM535 or HM536 (detection)
Preservatives	0.01% Sodium Azide