

Product Datasheet

Procalcitonin (PCT) Antigen, Recombinant ≥95% LA368

Product Name	Procalcitonin (PCT) Antigen, Recombinant ≥95%
Catalog Number	LA368
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (≤ 1 week): 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Purity	≥95% by SDS-PAGE
Buffer	50 mM Tris-HCl, pH 8.0 with 150 mM NaCl and 10% glycerol
Molecular Weight	13 kDa
Matched Pairs	HM563, HM564, HM698 or HM699
Expression System	E. coli
Preservatives	0.05% ProClin 300