

Product Datasheet

Troponin I, Cardiac Antigen, Recombinant >95% LA369L

Product Name Troponin I, Cardiac Antigen, Recombinant >95%

Catalog Number LA369L

Tested Applications ELISA, WB

Brand EastCoast Bio

Form/Appearance Lyophilized Purified Recombinant Antigen (His-tag)

Storage Short Term (< 1 week): 2-8°C. Long Term: ≤-20°C. After reconstitution: -20°C. Avoid repeated freezing and thawing.

Note Reconstitute with water to 4 mg/ml.

Purity >95% by SDS-PAGE

Molecular Weight 26 kDa

Expression System E. coli