

## Product Datasheet

### Fatty Acid Binding Protein (H-FABP) Antigen, Recombinant >95% LA370

**Product Name** Fatty Acid Binding Protein (H-FABP) Antigen, Recombinant >95%

**Catalog Number** LA370

**Tested Applications** ELISA, WB

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (His-tag)

**Storage** Short Term ( $\leq 1$  week): 2-8°C. Long Term:  $\leq -20^\circ\text{C}$ . Avoid repeated freezing and thawing.

**Note** Amino Acid Sequence:  
MVDAFLGTWKLVDKFNFDYMKSLGVGFATRQVASMTKPTTIIKNGDILTLKTHSTFKNTEISFKLGVFEDETTADDRK

**Purity** >95% by SDS-PAGE

**Buffer** 135 mM NaCl, 2.7 mM KCl, 1.5 mM  $\text{KH}_2\text{PO}_4$ , 8 mM  $\text{K}_2\text{HPO}_4$ , pH 7.2

**Molecular Weight** 15.9 kDa

**Expression System** E. coli