

Product Datasheet

Myoglobin (Mb) Antigen, Recombinant ≥95% LA372

Product Name Myoglobin (Mb) Antigen, Recombinant ≥95%

Catalog Number LA372

Tested Applications LF, ELISA, CLIA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage ≤-20°C. Avoid repeated freezing and thawing.

Purity ≥95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Matched Pairs HM687 or HM688

Expression System E. coli

Preservatives 0.09% Sodium Azide