

## Product Datasheet

### **Creatine Kinase BB (CK-BB) Antigen, Recombinant >95% LA374**

**Product Name** Creatine Kinase BB (CK-BB) Antigen, Recombinant >95%

**Catalog Number** LA374

**Tested Applications** ELISA, WB

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (His-tag)

**Storage** Short Term ( $\leq 1$  week): 2-8°C. Long Term:  $\leq -20^\circ\text{C}$ . Avoid repeated freezing and thawing.

**Note** Amino Acid Sequence:  
MPFSNSHNALKLRFPAEDEFPDLSAHNNHMAKVLTPELYAELRAKSTPSGFTLDDVIQTGVDNPGHPYIMTVGCVAGDE

**Purity** >95% by SDS-PAGE

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.4

**Molecular Weight** 44 kDa

**Expression System** E. coli