

## Product Datasheet

### Ebola Virus Antigen, Recombinant >95% LA382

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Ebola Virus Antigen, Recombinant >95%                      |
| <b>Catalog Number</b>      | LA382  |
| <b>Tested Applications</b> | LF, ELISA  |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Recombinant Antigen  |
| <b>Storage</b>             | -20°C. Avoid repeated freezing and thawing.                |
| <b>Purity</b>              | >95% by SDS-PAGE   |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Arginin |
| <b>Specificity</b>         | Ebola Virus (Zaire) IgG and IgM                            |
| <b>Expression System</b>   | E. coli  |
| <b>Preservatives</b>       | 0.02% Sodium Azide   |