

## Product Datasheet

### Ebola Virus Antigen, Recombinant >95% LA382

<b>Product Name</b>	Ebola Virus Antigen, Recombinant >95%
<b>Catalog Number</b>	LA382
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Recombinant Antigen
<b>Storage</b>	-20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Arginin
<b>Specificity</b>	Ebola Virus (Zaire) IgG and IgM
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.02% Sodium Azide