
Product Datasheet

Pinaca-BSA Antigen LA384

Product Name	Pinaca-BSA Antigen
Catalog Number	LA384
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	20 mM Phosphate Buffered Saline, pH 7.5
Matched Pairs	HM594 (detection)
Preservatives	0.01% Sodium Azide