

## Product Datasheet

### UR-144-BSA Antigen LA385

**Product Name** UR-144-BSA Antigen

**Catalog Number** LA385

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Liquid; Antigen-BSA Conjugate

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Buffer** 20 mM Phosphate Buffered Saline, pH 7.5

**Matched Pairs** HM595 or HM596 (detection)

**Preservatives** 0.01% Sodium Azide