

Product Datasheet

Anti-Mullerian Hormone (AMH) Antigen, Recombinant >90% LA387

Product Name Anti-Mullerian Hormone (AMH) Antigen, Recombinant >90%

Catalog Number LA387

Tested Applications LF, ELISA, CLIA

Brand EastCoast Bio

Form/Appearance Recombinant Antigen (His-tag)

Storage -20°C. Avoid repeated freezing and thawing.

Purity ≥95%

Buffer 50 mM Tris, pH 8.0 with 0.15M NaCl and 10% glycerol

Molecular Weight 32 kDa

Matched Pairs HM605, HM606, HM618, HM619, HM620 or HM621

Expression System E. coli

Preservatives 0.05% Proclin 300