

Product Datasheet

Anti-Mullerian Hormone (AMH) Antigen, Recombinant >90% LA388

Product Name	Anti-Mullerian Hormone (AMH) Antigen, Recombinant >90%
Catalog Number	LA388
Tested Applications	ELISA, CLIA
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen (full length with His-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	50 mM Tris, pH 8.0
Matched Pairs	HM607 or HM608
Expression System	E. coli
Preservatives	None