

Product Datasheet

Norovirus genogroup I (P domain) Antigen, Recombinant >95% LA391

Product Name	Norovirus genogroup I (P domain) Antigen, Recombinant >95%
Catalog Number	LA391
Tested Applications	ELISA, Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (P domain 225-520 with His-tag at N-terminus)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline with 25 mM arginine
Molecular Weight	30 kDa
Matched Pairs	HM637 or HM884
Expression System	E. coli
Preservatives	0.05% Sodium Azide