

Product Datasheet

Norovirus genogroup II (P domain) Antigen, Recombinant 95% LA392

Product Name	Norovirus genogroup II (P domain) Antigen, Recombinant 95%
Catalog Number	LA392
Tested Applications	ELISA, Control
Immunogen	Norovirus GII.4 VLP
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (P domain 225-520 with His-tag at N-terminus)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	PBS with 25 mM Potassium Carbonate, pH 9.0
Specificity	Norovirus GII
Molecular Weight	30 kDa
Matched Pairs	HM638 or HM885
Expression System	E. coli
Preservatives	0.05% Sodium Azide