

Product Datasheet

Helicobacter pylori (OMP) Antigen, Recombinant ≥90% LA396

Product Name Helicobacter pylori (OMP) Antigen, Recombinant ≥90%

Catalog Number LA396

Species/Host E. coli

Tested Applications LF, ELISA, TIA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage 2-8°C

Purity ≥90%

Buffer Phosphate buffer with additives and salt, pH 7.0

Matched Pairs HM629

Expression System E. coli

Preservatives 0.09% Sodium Azide