

Product Datasheet

Clostridium difficile Glutamate Dehydrogenase (GDH) Antigen, Recombinant >95% LA398

Product Name	Clostridium difficile Glutamate Dehydrogenase (GDH) Antigen, Recombinant >95%
Catalog Number	LA398
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	2-8°C
Purity	>95%
Buffer	Tris Buffer, pH 7.6
Molecular Weight	51 kDa
Matched Pairs	HM631
Expression System	E. coli
Preservatives	0.09% Sodium Azide