

## Product Datasheet

### SS-A Antigen, Recombinant LA412

<b>Product Name</b>	SS-A Antigen, Recombinant
<b>Catalog Number</b>	LA412
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (52 kD)
<b>Storage</b>	≤ -70°C
<b>Buffer</b>	20 mM Phosphate Buffered Saline, pH 6.0 with 500 mM NaCl and 3M Urea
<b>Expression System</b>	E. coli