

---

## Product Datasheet

### SCL-70 Antigen, Recombinant LA416

<b>Product Name</b>	SCL-70 Antigen, Recombinant
<b>Catalog Number</b>	LA416
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	≤-70°C
<b>Buffer</b>	20 mM Phosphate Buffer, pH 7.8 with 500 mM NaCl
<b>Expression System</b>	E. coli