
Product Datasheet

Escherichia coli (E. coli) O157 Antigen, Native LA420

Product Name Escherichia coli (E. coli) O157 Antigen, Native

Catalog Number LA420

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Native Inactivated Antigen

Storage Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Matched Pairs HM649

Preservatives 0.09% Sodium Azide