

## Product Datasheet

### Dengue Virus (DENV-Env Serotype 2) Antigen, Recombinant >90% LA421

<b>Product Name</b>	Dengue Virus (DENV-Env Serotype 2) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA421
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	-20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>90%
<b>Buffer</b>	20 mM Phosphate Buffered Saline
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.01% ProClin 300