

## Product Datasheet

### Salmonella typhi Antigen, Native LA424

**Product Name** Salmonella typhi Antigen, Native

**Catalog Number** LA424

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Liquid; Inactivated Native Antigen

**Storage** Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Source** S. typhi culture

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.4

**Matched Pairs** HM174, HM037, HM038, HM039, or HM040

**Preservatives** 0.09% Sodium Azide