
Product Datasheet

Salmonella paratyphi B Antigen, Native LA426

Product Name	Salmonella paratyphi B Antigen, Native
Catalog Number	LA426
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Inactivated Native Antigen
Storage	Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	S. paratyphi B culture
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM037, HM038, HM039, or HM040
Preservatives	0.09% Sodium Azide