

---

## Product Datasheet

### Corticosterone-3-BSA Antigen LA431

<b>Product Name</b>	Corticosterone-3-BSA Antigen
<b>Catalog Number</b>	LA431
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	HM655 or HM656 (detection)
<b>Preservatives</b>	0.05% Sodium Azide