

Product Datasheet

Lysergic acid diethylamide (LSD)- BSA Antigen LA434

Product Name	Lysergic acid diethylamide (LSD)- BSA Antigen
Catalog Number	LA434
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	0.015M Phosphate, pH 7.2 with 0.15M NaCl
Preservatives	0.05% Sodium Azide