

Product Datasheet

Lysergic acid diethylamide (LSD)- BSA Antigen LA434

| | |
|----------------------------|---|
| Product Name | Lysergic acid diethylamide (LSD)- BSA Antigen |
| Catalog Number | LA434 |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid; Antigen-BSA conjugate |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Buffer | 0.015M Phosphate, pH 7.2 with 0.15M NaCl |
| Preservatives | 0.05% Sodium Azide |