

## Product Datasheet

### Hepatitis B Core Antigen (HBcAg) Antigen, Recombinant >90% LA439

**Product Name** Hepatitis B Core Antigen (HBcAg) Antigen, Recombinant >90%

**Catalog Number** LA439

**Tested Applications** LF, ELISA, WB

**Brand** EastCoast Bio

**Form/Appearance** Recombinant Antigen

**Storage** Short Term: 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.

**Purity** >90%

**Buffer** 20 mM Tris-HCl, pH 8.0 with 1 mM EDTA

**Molecular Weight** 15 kDa

**Expression System** E. coli

**Preservatives** 0.1% Proclin 300