

## Product Datasheet

### Alpha-fetoprotein (AFP) Antigen, Recombinant >95% LA443

**Product Name** Alpha-fetoprotein (AFP) Antigen, Recombinant >95%

**Catalog Number** LA443

**Tested Applications** Control

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen

**Storage** ≤-20°C. Avoid repeated freezing and thawing.

**Purity** >95%

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.4

**Molecular Weight** 68 kDa

**Expression System** Mammalian Cell

**Preservatives** 0.05% Proclin 300