

Product Datasheet

Carcinoembryonic Antigen (CEA), Recombinant >95% LA444

Product Name Carcinoembryonic Antigen (CEA), Recombinant >95%

Catalog Number LA444

Tested Applications Control

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage ≤-20°C. Avoid repeated freezing and thawing.

Purity >95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Molecular Weight 100-180 kDa

Expression System Mammalian Cell

Preservatives 0.05% Proclin 300