

Product Datasheet

Creatine Kinase MB (CK-MB) Antigen, Recombinant >95% LA445

Product Name Creatine Kinase MB (CK-MB) Antigen, Recombinant >95%

Catalog Number LA445

Tested Applications Control

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage ≤-20°C. Avoid repeated freezing and thawing.

Purity >95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Molecular Weight 88 kDa

Expression System E. coli

Preservatives 0.05% Proclin 300