

Product Datasheet

ST2 Antigen, Recombinant >95% LA446

Product Name	ST2 Antigen, Recombinant >95%
Catalog Number	LA446
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	≤-20°C. Avoid repeated freezing and thawing.
Purity	≥95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM696 or HM697
Expression System	Mammalian Cell
Preservatives	0.02% Proclin 300