

## Product Datasheet

### Lipoprotein a (Lp-a) Antigen, Native >90% LA447

<b>Product Name</b>	Lipoprotein a (Lp-a) Antigen, Native >90%
<b>Catalog Number</b>	LA447
<b>Tested Applications</b>	Immunoassays, Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Native Antigen
<b>Storage</b>	≤-20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90% by SDS-PAGE
<b>Source</b>	Human Serum
<b>Buffer</b>	0.02M PBS, pH 7.4
<b>Preservatives</b>	0.05% ProClin 300