

Product Datasheet

Gastrin 17 (G-17) Antigen, Synthetic >95% LA451L

Product Name	Gastrin 17 (G-17) Antigen, Synthetic >95%
Catalog Number	LA451L
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Lyophilized Purified Synthetic Antigen
Storage	≤ -20°C. Avoid repeated freezing and thawing.
Note	Reconstitute with 10 mM PBS, pH 7.4
Purity	>95% by SDS-PAGE
Source	Synthetic
Matched Pairs	HM774, HM775, HM776