

## Product Datasheet

### Gastrin 17 (G-17) Antigen, Synthetic >95% LA451L

<b>Product Name</b>	Gastrin 17 (G-17) Antigen, Synthetic >95%
<b>Catalog Number</b>	LA451L
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized Purified Synthetic Antigen
<b>Storage</b>	≤ -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Reconstitute with 10 mM PBS, pH 7.4
<b>Purity</b>	>95% by SDS-PAGE
<b>Source</b>	Synthetic
<b>Matched Pairs</b>	HM774, HM775, HM776