

Product Datasheet

Malaria (pv-LDH) Antigen, Recombinant >90% LA454

Product Name	Malaria (pv-LDH) Antigen, Recombinant >90%
Catalog Number	LA454
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	20 mM Tris-HCl, pH 8.0 with 0.1 mM EDTA
Molecular Weight	40 kDa
Expression System	E. coli