

Product Datasheet

Acrosomal Protein SP-10 (ACRV1) Antigen, Recombinant >90% LA457

Product Name Acrosomal Protein SP-10 (ACRV1) Antigen, Recombinant >90%

Catalog Number LA457

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage $\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Matched Pairs HM1185-HM1192

Expression System Mammalian Cell

Preservatives 0.05% Sodium Azide