
Product Datasheet

Herpes Simplex Virus 2 (HSV 2-gB) Antigen, Recombinant >90% LA465

Product Name Herpes Simplex Virus 2 (HSV 2-gB) Antigen, Recombinant >90%

Catalog Number LA465

Tested Applications ELISA, WB

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage Short Term (<1 week): 2-8°C. Long Term: -20°C

Purity >90% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.6

Molecular Weight 80 kDa

Expression System E. coli